

REACH Perimed declaration regarding (EC) 1907/2006

Perimed AB declares that to the best of our knowledge, including information from our suppliers, our products comply with the European Regulation (EC) 1907/2006 requirements regarding Registration, Evaluation, Authorization and Restriction of Chemicals (REACH).

Compliance with REACH requires that:

- Manufacturers and importers must register chemical substances that they produce or import in quantities over 1 t/year
- Manufacturers must inform customers if their products contains substances from the SVHC Candidate List in a concentration above 0,1% weight by weight (w/w)
- Substances in Annex XVII are not allowed to be used if the conditions of restrictions apply to Perimeds products

Perimed strives to substitute all Substances of Very High Concern (SVHC) that are found in our products in a concentration above 0.1% (w/w), when technically possible.

Based on our investigations and information from our suppliers we have identified the following substances from the SVHC candidate list, in a concentration above 0.1% (w/w), in our products/articles:

PeriCam PSI NR/HR (90-00070 / 90-00071)

Articles affected, Product name	Articles affected, Part number	Affected part, name & number	Substance name	CAS no. (where available)	RoHS (2011/65/EU) exemption (if applicable) Specified only for products CE-marked by Perimed
High-end Laptop PC for PeriCam PSI	93-00311	Laptop Dell Precision MW 7550 (79-0000270-01)	Several	See 79-0000270- 01 Attachment 4 REACH	N/A

Articles affected, Product name	Articles affected, Part number	Affected part, name & number	Substance name	CAS no. (where available)	RoHS (2011/65/EU) exemption (if applicable) Specified only for products CE-marked by Perimed
High-end Laptop PC for PeriCam PSI	93-00311	Laptop Dell Precision MW 7560 (79-0000270-02)	Several	See 79-0000270- 02 Attachment 7 REACH	N/A
High-end Desktop PC for PeriCam PSI	93-00310	Desktop Dell (79-0000280)	Several	See 79-0000280 Attachment 4 REACH	N/A
High-end Panel PC for PeriCam PSI	93-00312	Panel PC HP (79-0000272)	Several	See 79-0000272 Attachment 6 REACH	N/A
Panel PC for PeriCam PSI	93-00309	Panel PC win 10 (79-0000223)	Several	See PeriFlux 6000 79-0000223	N/A
Isolation Station	93-00313	Isolation Station (79-0000176)	Lead	7439-92-1	Annex III; 6c, 7a
PeriCam PSI NR/HR	90-00070 90-00071	PCBA PC1420 (41-00683)	Lead	7439-92-1	Annex III; 6c, 7a
PeriCam PSI NR/HR	90-00070 90-00071	PCBA PC1377 (41-00245)	Lead	7439-92-1	Annex III; 7a
Vexio-cart Plus for Imager	93-00307	Power cable US/Canada 15 Feet	Lead	7439-92-1	Annex III; 6c, 8b
Vexio-cart Plus for Imager, height Adjustable	93-00308	(73-0600040)			

**PeriFlux 6000 Systems (93-00200, 93-00202, 93-00205), and
PeriFlux 6000 tcpO2 stand-alone (93-00204)**

Articles affected, Product name	Articles affected, Part number	Affected part, name (part number)	Substance name	CAS no. (where available)	RoHS (2011/65/EU) exemption (if applicable) Specified only for products CE-marked by Perimed
Alt. 1 Panel PC *	93-00179*	HP Panel PC win 10 (79-0000223)	Several	See 79-0000223 Attachment 3 REACH	N/A
Alt. 2 Panel PC *	93-00179*	Fujitsu Panel PC win 10 (79-0000281)	Several	See 79-0000281 Attachment 2 REACH	N/A
PF 6001 Main Unit	90-00056 90-00067	PCBA PC1384 (41-00373-01)	1: Lead 2: Cadmium 3: 1, 2 dimethoxyethane; ethylene glycol dimethyl ether (EGDME)	1: 7439-92-1 2: 7440-43-9 3: 110-71-4	Annex III; 6c, 7a, 7c-I, 8a
PF 6001 Main Unit	90-00056 90-00067	Center module (41-00352-01) including COMe-mBTi10 E3825 2GB, PERIMED (79-0000292)	1.Frits, Chemicals (Pb containing) 2. Lead 3. Lead monoxide (lead oxide) 4. Diboron trioxide 5.4,4'Isopropylidenediphenol (BPA)	1: 65997-18-4 2: 7439-92-1 3: 1317-36-8 4: 1303-86-2 5: 80-05-7	Annex III; 7a & 7c-I
PF 6001 Main Unit	90-00056 90-00067	Center module (41-00352-01) including COMe-mEL10 x6212RE 4GB, PERIMED (79-0000293)	1.Lead 2.Frits Chemicals (Pb containing) 3.Lead titanium trioxide	1: 7439-92-1 2: 65997-18-4 3: 12060-00-3	Annex III; 7a & 7c-I

Articles affected, Product name	Articles affected, Part number	Affected part, name (part number)	Substance name	CAS no. (where available)	RoHS (2011/65/EU) exemption (if applicable) Specified only for products CE-marked by Perimed
PF 6010 LDPM/Temp Unit	90-00059	PCBA PC1412 (41-00471)	Lead	7439-92-1	Annex III; 7a & 7c-I
PF 6010 LDPM/Temp Unit	90-00059	ST-connector, ceramic (77-0001401)	Lead	7439-92-1	Annex III; 6c
PF 6010 LDPM/Temp Unit	90-00059	Glue EPO-TEK 302-3M (87-1000056)	4-Nonyl Phenol	84852-15-3	N/A
PF 6040 tcpO2/pCO2 Unit	90-00057	PCBA PC1382 (41-00362-02)	Lead	7439-92-1	Annex III; 7a & 7c-I
PF 6050 Pressure Unit	90-00060	80-0000249 Air pressure connector	Lead	7439-92-1	Annex III; 6c
PF 6050 Pressure Unit	90-00060	80-0000250 Air pressure connector	Lead	7439-92-1	Annex III; 6c
PF 6050 Pressure Unit	90-00060	PCBA PC1415 (41-00468)	Lead	7439-92-1	Annex III; 7a & 7c-I
PeriFlux 6000 Systems*	*	Adapter Cord (79-0000189)	Lead	7439-92-1	Annex III; 6c
PeriFlux 6000 Systems* Isolation Transformer with Network Isolator Isolation Transformer	* 93-00209 93-00210	Isolation Station (79- 0000176)	Lead	7439-92-1	Annex III; 6c, 7a
PeriFlux 6000 Systems*	*	Power cable US/Canada 15 Feet (73-0600040)	Lead	7439-92-1	Annex III; 6c, 8b

*Included in PeriFlux 6000 Systems (93-00200, 93-00202, 93-00205)

Perilont (93-00123, 93-00124)

Articles affected, Product name	Articles affected, Part number	Affected part, name (part number)	Substance name	CAS no. (where available)	RoHS (2011/65/EU) exemption (if applicable)
PF 751 Perilont Power Supply	41-00321	PCBA PC1381 Perilont (41-00323)	Lead	7439-92-1	Annex III; 7a & 7c-I

Probes

Articles affected, Product name	Articles affected, Part number	Affected part, name & number	Substance name	CAS no. (where available)	RoHS (2011/65/EU) exemption (if applicable)
PROBE 418 Master Probe PROBE 404-1 Suturable Angled Probe PROBE 407 Small Straight Probe PROBE 481 Iontophoresis Probe for LDPM PROBE 404-NI Non-invasive Angled Probe PROBE 457 Small Angled Thermostatic Probe PROBE 411, Needle Probe microm	91-00053 91-00219 91-00023 91-00001 91-00232 91-00015 91-00050	Y-Joint cable branching (13-00879-01) including PVC Black (87-1000077)	dioctyl-tin-combinations	15571-58-1	N/A


Articles affected, Product name	Articles affected, Part number	Affected part, name & number	Substance name	CAS no. (where available)	RoHS (2011/65/EU) exemption (if applicable)
PROBE 418 Master Probe PROBE 404-1 Suturable Angled Probe PROBE 407 Small Straight Probe PROBE 481 Iontophoresis Probe for LDPM PROBE 404-NI Non-invasive Angled Probe PROBE 457 Small Angled Thermostatic Probe PROBE 411, Needle Probe microm	91-00053 91-00219 91-00023 91-00001 91-00232 91-00015 91-00050	ST-connector, ceramic (77-0001401)	Lead	7439-92-1	Annex III; 6c
PROBE 418 Master Probe PROBE 404-1 Suturable Angled Probe PROBE 407 Small Straight Probe PROBE 481 Iontophoresis Probe for LDPM PROBE 404-NI Non-invasive Angled Probe PROBE 457 Small Angled Thermostatic Probe PROBE 411, Needle Probe microm	91-00053 91-00219 91-00023 91-00001 91-00232 91-00015 91-00050	Glue EPO-TEK302-3M (87-1000056)	4-Nonyl Phenol	84852-15-3	N/A

The above information is to the best of our knowledge accurate. Perimed will continue to monitor our supply chain for continued compliance with REACH.

Perimed AB

Järfälla, DATE

2022-03-16



Lena Åredal
President